

First Toric ARTIFLEX implanted in USA



Left to right: Aleksandra Stopinska, C.O.T., Andrew C. Shatz, M.D., Jean Pak, OPHTEC USA, Victoria Groeniger, R.N.



We are proud to announce that Dr Andrew Shatz, SightTrust Eye Institute, Sunrise, Florida, implanted the first Ophtec ARTIFLEX Toric Lens in the United States on August 22, 2013. Dr Shatz was granted special approval by the FDA to use this amazing 'implantable contact lens' to help a patient who had become unable to tolerate contact lenses or glasses and was legally blind without them.

Prior to surgery, the patient was unable to see things that were more than a few inches from the face without strong prescription lenses. The patient's vision at the 1-day post-operative visit was 20/25 at distance and 20/20 at near. This patient was thrilled and will now be able to pursue a lifelong dream of serving in the United States Air Force!

ARTIFLEX Toric IOL



- ✓ IRIS FIXATION
- ✓ SMALL INCISION
- ✓ NO ROTATION POST OP
- ✓ LARGE OPTICAL ZONE

Optic:	6.0 mm
Haptics:	Iris Fix®
Overall diameter:	8.5 mm
Dioptric range:	-1.0 D to -13.5 D in combination with a cylinder of -1.0 D to -5.0 D

More Ophtec (P)IOLs through www.ophtec.com

Find Ophtec in your region:

The Netherlands (head office)

Ophtec BV
Schweitzerlaan 15
9728 NR Groningen,
The Netherlands
Tel: +3150 5251944 | Fax +31505254386
info@ophtec.com

Asia Pacific

Ophtec Asia-Pacific CO, Ltd.
Flat 9D Willow Mansion 22
Taikoo wan Road Taikoo Shing,
Hong Kong
Tel: +852 28871762 | alex@ophtec-ap.com

Germany

Ophtec GmbH
Drögensee 12 a, 22397 Hamburg,
Germany
Tel: +49 40 60096 978
a.lindstaedt@de.ophtec.com

Japan

Ophtec Japan Inc.
2-7-29, Kitaoyama, Minato-ku,
Tokyo 107-0061
Tel: +81 35919 4366 | yamada@ophtec-ap.com

North America

Ophtec USA Inc.
6421 Congress Ave., Suite 112, Boca Raton
FL 33487
USA
Tel: +1 561 989-8767
rick.mccarley@usa.ophtec.com

Portugal

Ophtec Portugal
Edificio Central Park
Rua Alexandre Herculano, Torre 1, 4 B
2795-240 Linda-a-Velha
Portugal
Tel: +351 215 900 801 | Fax: +351 210 936 453
s.bayan@prt.ophtec.com

South Africa

Ophtec South Africa (Pty) Ltd
Visiomed Office Park, Building 1 Unit 4,
269 Beyers Naude Drive, Blackheath Ext.1
Johannesburg, South Africa
Tel: +27 105903135 | deon@ophtec.co.za

South Korea Co, Ltd

Ophtec South Korea
7th FL. YungJeon Building 154-10,
SamSeoung-Dong, GangNam-Gu, Seoul, 135-879
South Korea
Tel: +82-2-508-0522 | y.m.goo@kr.ophtec.com

Spain

Lio Ophtec España SL
C/ Azalea, 1 Edif. B Planta 1
Oficina 3 Miniparc I-El Soto
28109 Alcobendas, Madrid
España
Tel: 900 993 174
y.hernandez@es.ophtec.com

